911 CANYON

880.01

Directive 93/42/EEC

EN

WWW.KONG.IT

Laden Sie die Übersetzung in Ihrer Sprache herunter - Download the translation in your language - Bájate la traducción en tu idioma - Télécharger la traduction dans votre langue - Scarica la traduzione nella tua lingua - 下载您语言版本的手册

TABLE OF CONTENTS

1 - SYMBOLS AND SUPPORT	4
1.1 Symbols	4
1.2 Technical assistance	4
2 - GENERAL INFORMATION	5
3 - TECHNICAL FEATURES	6
3.1 Terminology and materials of the parts	6
3.2 Dimensions	7
3.3 Range	7
3.4 Optional components and spare parts	10
4 - SPECIFIC INFORMATION	11
4.1 Intended use	11
4.2 Spoke Pattern	11
4.3 Positioning and immobilising the patient	14
4.4 Transporting the patient	17
4.5 Transporting of the stretcher	19
5 - MAINTENANCE AND REPAIR	20
5.1 General Information	20
5.2 Maintenance	20
5.3 Repair	20
6 - STORAGE	21
7 - CHECKS, INSPECTIONS AND SERVICING	22
7.1 Checks	22
7.2 Inspections	22
7.3 Revisions	22
8 - DURATION AND WARRANTY	23
8.1 Duration	23
8.2 Disposal	23
8.3 Warranty	23
8.4 Law obligations	23
9 - LABELLING AND SYMBOLS	24
9.1 Labelling	24
9.2 Symbols	24
10 - DOCUMENTS	25
10.1 Inspections and servicing register	25
10.2 Maintenance and repair records	26
10.3 Declaration of conformity (facsimile)	27
11 - LIST OF APPLIED AND REFERENCE STANDARDS	28
11.1 Applied standards	28
11.2 Standards used as reference	28

CHAPTER

1

SYMBOLS AND SUPPORT

1.1 SYMBOLS

For the sake of comfort and clarity while reading this manual, below are the symbols used for the management of important warnings to ensure correct and safe use of the device.

REQUIREMENT FOR PROPER USE

Identifies the presence of information for proper use of the device.

INFORMATION REQUIREMENT

Identifies the presence of useful and general information which guides the reader towards conscientious use of the device and/or performance of actions.

Shows that the product is made, designed and manufactured in compliance with the safety requirements (RES) laid down by the Medical Devices Directive 93/42/CEE (Class I medical devices, in compliance with the classification regulation 1 as shown in annex IX).

1.2 SUPPORT

For information please contact Kong Customer Support Service by:

- telephone 0039 0341 630506,
- fax 0039 0341 641550,
- email: safetycare@kong.it,

or write to KONG S.p.A. – Via XXV Aprile, 4 – 23804 Monte Marenzo LC - ITALY.

To facilitate assistance operations, please always communicate or state the serial number (SN) indicated on the label fixed to the Medical Device.

CHAPTER 2

2 GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the device. The information is about the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers advice on use, this information shall not be treated as a user manual under actual conditions of use.

WARNINGS AND LIMITATIONS OF USE:

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, dropping, wear, chemical contamination, exposing textile/plastic components/

devices to temperature below -30°C or above $+50^{\circ}\text{C}$ and metal components/ devices to temperatures exceeding 100°C are some of the examples of causes that can reduce, limit and end the life of the device,

- prior to any rescue operation, ensure not to exceed the capacity indicated in paragraph 3.3,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or accessories different from the ones indicated in paragraph 3.4 can be,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase risks related to cutting and abrasion on textile and synthetic devices.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device in the natural environment may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/expert/qualified people.

Warning: not suitable for use in an ATEX \ environment (Directive 94/9/EC)

CHAPTER 3

TECHNICAL FEATURES

3.1 TERMINOLOGY AND MATERIALS OF THE PARTS

- A – Aluminium alloy frame, can be disassembled into two parts,
- B – Protective shells made of GRP,
- C – Aluminium alloy handles,
- D – Polyurethane padding,
- E – Waterproof canvas in trilaminare,
- F – Aeration valve in ABS,
- G – Cranial protection in Lexan®,
- H – Polyester shoulder straps,
 - Mattresses in polyurethane foam,
 - Polyester webbing.


Fig. 1

3.2 DIMENSIONS

Length without handles:	220 cm
Length with handles:	350 cm
Width:	60 cm
Overall weight (fig.1):	23 kg

3.3 CARRYING CAPACITY

The table below defines the capacities of the “911 CANYON” stretcher, based on the method of use deemed most appropriate by the rescue operator, the outcome of the laboratory tests and recommended safety coefficients.

How to use	Breaking load	Safety factor	Capacity
Carriage with handles/sun loops on the side of the frame (fig. 2-3)	> 450 kg	1:3	150 kg
Lifting and lowering using rescue helicopters (fig. 4) <i>by suspending</i>	> 1500 kg	1:10	150 kg
Transport by cable car and airlines (figg. 4-5) <i>by suspending</i>	> 1500 kg	1:10	150 kg
Vertical descents and recoveries (fig. 6) <i>by rope</i>	> 900 kg	1:6	150 kg

Prior to any rescue operation, ensure not to exceed the indicated capacity.


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6

3.4 OPTIONAL COMPONENTS AND SPARE PARTS

3.4.1 Optional components

The “911 CANYON” stretcher is complete with all the optional components for the use described in this manual.

3.4.2 Spare parts

- D - Paddings for handles,
- G - Cranial protection,
- H - Shoulder straps for stretcher,
- I - Suspended,
- L - Rope bag.


Fig. 7

CHAPTER 4

SPECIFIC INFORMATION

4.1 INTENDED USE

The “**911 CANYON**” stretcher is a medical device particularly suitable for recovering and transporting a patient in a gorge, even if immobilized on the “**VACUUM**” vacuum mattresses. Decisions concerning moving and immobilising the patient, as well as the duration, method to be used and combination with other devices, shall be taken and executed by expert and trained personnel only.

4.2 ASSEMBLY

Assemble the two parts of the structure by means of the following operations:

- rotate the two safety rods (Fig. 8),
- lift and turn the two locking rods (Fig. 9),

Fig. 8


Fig. 9

c) connect the two parts of the structure by means of the frame pins,


Fig. 10

- d) reposition the two locking rods in their respective seats (Fig. 11),
- e) reposition the two safety rods (fig. 12) and check the tightness of the coupling,
- f) Screw the two rings together completely (fig. 13).

Fig. 11


Fig. 12



Fig. 13

g) attach the watertight sheet to the lower part of the structure (feet) by means of the automatic buckle.


Fig. 14

h) lay down the watertight sheet.


Fig. 15

4.3 POSITIONING AND IMMOBILISING THE PATIENT

With the scoop stretcher positioned on a flat and stable surface:

- a) release the buckles and open the outer straps,
- b) open the sheet completely, holding both edges of the zipper and avoiding tearing,


Fig. 16

- c) connect the cranial protection to the structure:
 - insert the fins into the slots provided (fig. 17),
 - close the nylon buckles and tension the webbing (fig. 18).

All fins must rest on the Lexan® worktop


Fig. 17


Fig. 18

d) release the buckles, open the locking harness, position and adjust it,


Fig. 19

e) place the patient in the stretcher,


Fig. 20

f) secure the patient with the harness,

- This procedure should also be applied in the case of patients immobilised on the "VACUUM" vacuum mattress,
- the harness can also be used in part,
- the patient's arms are held by the tensioning of the outer straps.

Fig. 21

g) close the watertight sheet with the hinges, attach the buckles of the outer straps and tighten them.


Fig. 22

4.4 TRANSPORTING THE PATIENT

The “911 CANYON” stretcher is suitable for lifting and transporting the patient using the handles (fig. 2), the side handles (fig. 3), suspension (fig. 4-5) and lifting/calculating the patient vertically using the static rope (fig. 6).

4.4.1 Connection with provided handles

a) Connect the handles to the stretcher frame, by using the coupling numbering,


Fig. 23

b) lock the connection by inserting the shaped plug in the hole.


Fig.24

4.4.2 Connection with suspension

- Connect the suspension to the side holes of the stretcher using the four connectors,
- screw the connector nuts in completely.



Fig. 25

The eyelets of the suspension allow you to adjust the position of the stretcher.

4.4.3 Connection with ropes

- Connect the ropes to the attachment points (quick links) at the ends of the stretcher using the connectors,
- screw in the connector nuts completely.


Fig. 26

4.5 TRANSPORTATION OF THE STRETCHER

Attach a shoulder strap to each part of the structure as defined in Figures 27-28 and 29.


Fig. 27


Fig. 28


Fig. 29

CHAPTER 5

MAINTENANCE AND REPAIR

5.1 GENERAL INFORMATION

The “**911 CANYON**” stretcher is made of material which is highly resistant to wear and external agents. Despite this, the use conditions require maintenance operations and even repair in particular conditions.

Maintenance and repair operations must be recorded in a special register, an example of which is indicated in chapter 10.

5.2 MAINTENANCE

Below are maintenance operations to be carried out by the user:

- a) cleaning: after every use, wash with warm (max. 40° C) drinking water, possibly adding a neutral detergent (e.g. Marseille soap). Rinse and leave to dry in a shaded place, away from direct sources of heat,
- b) disinfection, when deemed necessary: soak the product in warm water in which sodium hypochlorite has been diluted at 1% for one hour, then rinse thoroughly using drinking water and leave to dry in a shaded place, away from direct sources of heat.

5.3 REPAIR

Repair shall be strictly carried out by the manufacturer.
The user is allowed to replace the parts mentioned in paragraph 3.4.2 only, with new and genuine parts.

CHAPTER 6

6 STORAGE

After cleaning, any disinfection and drying operations, keep the product and its accessories in a place that is dry (relative humidity 40-90%), cool (temperature 5-40° C) and safe (avoid UV radiation), chemically neutral (absolutely avoid salty and/or acid environments), away from sharp edges, sources of heat, humidity, corrosive substances or other possible jeopardising conditions.

Do not store this device wet!

CHAPTER 7

CHECKS, INSPECTIONS AND SERVICING

7.1 CHECKING

Before and after every use, check the device and ensure that:

- it is appropriate for the use it was intended for,
- it has not been subjected to mechanical deformations and it neither has cracks nor shows signs of wear,
- the textile parts do not have cuts, burns, chemical products residues, excessive hair, wear, in particular check areas at contact with the metal components (buckles, rings, etc.);
- the markings, including labels, are legible.

7.2 INSPECTION

The device must be inspected at least once a year, starting from the date of first use, by expert personnel trained and approved by KONG S.p.A. The date of the first use and the outcome of the inspection must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

7.3 CAPACITY

The device must be serviced by KONG S.p.A., or by expert personnel trained and approved by KONG S.p.A. if:

- any malfunctions are detected,
- the result of the pre and post use check or of the inspections is negative.

The outcome of the servicing must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

The revised device is guaranteed for one year from the revision date.

CHAPTER 8

DURATION AND WARRANTY

8.1 DURATION

Carefully read the “Warnings and limitations of use” paragraph in chapter 2. The duration of this device depends on the outcome of the at least one-yearly inspection (paragraph 7.2) and possible servicing. The suspension, even if not used, has a maximum duration of 10 years from the date of manufacture indicated on the label.

Important: any devices that have not passed the pre-use, post-use and periodic inspection must be eliminated and made unusable.

8.2 DISPOSAL

Follow the rules in force in the country of use and the waste disposal procedures of the reference hospital facility for proper disposal.

8.3 WARRANTY

The manufacturer guarantees the compliance of the device with the regulations in force at the time of production. Warranty covering defects is limited to defects regarding raw materials and manufacturing defects. The warranty does not cover normal wear, oxidation, damage caused by non-compliant use and/or in competition, improper maintenance, transportation, preservation or storage. Warranty will be immediately deemed null and void should the device be modified or tampered with. Validity corresponds to the legal validity of the country where the device was sold, as from the date of sale, by the manufacturer. No claims shall be raised against the manufacturer once the aforementioned period has expired. Any repair or replacement request under warranty must be accompanied by a proof of purchase. Should the defect be recognised, the manufacturer undertakes to repair or, replace or refund the device at his discretion. The manufacturer shall not be held liable beyond the invoice price of the device in any case whatsoever.

8.4 LAW OBLIGATIONS


Professional and free time activities are often regulated by special national laws which can impose limits and/or obligations concerning the use of these devices. The user has the obligation to know and apply such laws which could provide for limits different from those indicated in this information.

CHAPTER 9

LABELLING AND SYMBOLS

9.1 LABELLING

Labelling example.

Any changes to the positioning of the symbols does not change the content thereof.

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Consult the user manual

Warning: safety-related information, refer to the user manual

Compliant with Council Directive 93/42/EEC et seq. as amended, Class I DM

10.3 DECLARATION OF CONFORMITY (FACSIMILE)

This device comes with its own EC Declaration of Conformity with a signed original copy. In case of loss or should the packaging be damaged, you may request replacement at: safetycare@kong.it by submitting the unique serial number [SN] indicated on the label of the device.

KONG s.p.a.

via XXV Aprile, 4 | 23804 Monte Marengo (LC) - ITALY | tel +39 0341.630506 | info@kong.it

certified UNI EN ISO 9001

DECLARATION OF CONFORMITY

Manufacturer:

Company name: **KONG S.p.A.**
Registered office: **Via XXV Aprile, 4 - 23804 Monte Marengo (LC)**
Operational offices: **Via XXV Aprile, 4 - 23804 Monte Marengo (LC)**
VAT Code: **IT 00703180166**

hereby declares that the:

Medical device: **911 CANYON**
Class: **I according to attachment IX regulation 1**
REF: **880.01**
SN (serial number): _____
Date of manufacture: _____
registered in the Ministry of Health data bank under number **1448760**,

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complies with the essential requirements laid down by the Italian Legislative Decree n° 46/97 and subsequent amendments and additions transposing the Medical Devices Directive 93/42/EEC and subsequent amendments and additions.

Monte Marengo, _____

KONG S.p.A.
The Legal Representative
Mr. Marco Bonaiti

Revision form 0 dated 01/08/2016

CHAPTER 11

LIST OF APPLIED AND REFERENCE STANDARDS

11.1 APPLIED STANDARDS

- EN 1041:2008
Information provided by the medical devices manufacturer
- EN ISO 14971:2012
Medical devices - Risk management implementation for medical devices
- EN ISO 15223-1:2016
Medical devices - symbols to be used on the labels of the medical device, in the labelling and information that must be provided - Part 1: General requirements
- EN 62366:2008
Medical devices – Application of engineering for the exploitative characteristics of medical devices

11.2 STANDARDS USED AS REFERENCE

- EN 1865-1:2010+A1:2015
Patient transportation equipment used in ambulances - Part 1: General systems for stretchers and patient transportation equipment
- EN ISO 13485:2016
Medical devices: quality management systems; requirements for regulatory purposes
- EN 13718-1:2008
Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances

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